



Low Input RNAseq libraries - SMART-Seq v4 cDNA (1/2 reaction)

SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing (Clontech Cat.# from 634888 to 634894)

SeqAmp™ DNA Polymerase (Clontech Cat.#638504)

- Prepare a 10X Reaction Buffer by mixing 19 ul of 10X Lysis Buffer and 1 ul of RNase Inhibitor
- Dilute the samples in Reaction Buffer:

		1/4 reaction	___rxns
10X Reaction Buffer	0.5 ul	0.3 ul	_____
3'SMART-Seq CDS P. IIA	<u>1</u> ul	<u>0.5</u> ul	_____
	1.5 ul	0.8 ul	
Sample (Cells or RNA)	4.75 ul	2.4 ul	
Total volume	6.25 ul	3.2 ul	

hot lid (80°C)

72°C - 3 min
4°C - forever

- Make master mix with 10% extra rxn:

		1/4 reaction	
5X First-Strand Buffer	2 ul	1 ul	_____
SMART-Seq v4 Oligo	0.5 ul	0.3 ul	_____
RNase Inhibitor	0.25 ul	0.2 ul	_____
SMARTScribe enzyme	<u>1</u> ul	<u>0.5</u> ul	_____
Total volume per rxn	3.75 ul	2 ul	

add to 6.25 ul - Total 10ul

hot lid (80°C)

42°C - 90 min
70°C - 10 min
4°C - forever

SAFE STOP OVERNIGHT

- Make master mix with 10% extra rxn:

		1/4 reaction	
2X SeqAmp Buffer	12.5 ul	6.3 ul	_____
PCR Primer IIA	0.5 ul	0.3 ul	_____
SeqAmp Polymerase	0.5 ul	0.3 ul	_____
Nuclease Free Water	<u>1.5</u> ul	<u>0.8</u> ul	_____
Total volume per rxn	15 ul	7.7 ul	

add to 10 ul - Total 25ul

95°C - 1 min
<u>cycles</u> X _____
98°C - 10 sec
65°C - 30 sec
68°C - 3 min
72°C - 10 min
4°C - forever

SAFE STOP OVERNIGHT

RNA	Cells	cycles
5000 pg	500	9
500 pg	50	12
50 pg	5	16
5 pg	1	19

- Add 25 ul of Ampure XP beads (1.0x), wash, dry and elute in 17 ul of Elution Buffer. 1/4 reaction (13 ul)
- Quantify with BioA DNA HS chip - range 200bp to 9000bp.

